Ultravision™ Visual Field Clearing System

Ultravision™ 5mm Trocar


- DE Ultravision™ 5 mm Trokar
- FR Trocart 5 mm Ultravision™
- ES Trocar de 5 mm de Ultravision™
- IT Trocar Ultravision™ 5 mm
- NL Ultravision™ 5 mm trocar
- PT Trocarte de 5 mm Ultravision™
- SV Ultravision™ 5mm troakar
- EL Τροκάπ Ultravision™ 5 χιλιοστών
- DA Ultravision™ 5 mm trokar
- FI Ultravision™ 5mm -troakaari
- PL Trokar 5 mm Ultravision™
- LT Ultravision™ 5 mm troakaras
- SK 5mm Trokár Ultravision™
- JA Ultravision™ 5mm トロカール

Instructions, Gebrauchsanleitung, Consignes, Instrucciones, Istruzioni,
Instrucții, Instruções, Bruksanvisningar, Οδηγίες, Vejledninger,
Käyttöohje, Instrukcje, Instrukcijos, Návod, 取扱説明書


Please read the Ultravision User Manual and these Instructions for Use before using this product. Failure to do so may lead to serious injury to patient and operator.

Rx
ONLY


①


②


③


4

5

6

Important

This is not a complete set of instructions. Users must read the Ultravision User Manual before using this product. This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

Device description

The Ultravision™ 5mm Trocar is intended for use with the Ultravision™ Visual Field Clearing System. The pack contains a sterile, single use device consisting of an Ionwand™, trocar, main obturator (BLUE) and Ionwand obturator (RED). The trocar features an internal seal system that maintains pneumoperitoneum when instruments are inserted or withdrawn. The seal accommodates instruments up to 5mm in diameter.

The Ultravision 5mm Trocar is designed, tested and manufactured for single use only. Reuse or reesterilization of this product is prohibited to avoid contamination, patient infection and/or device malfunction and subsequent patient injury.

Contents (Figure 1)

1. Ionwand™
2. Trocar
3. Main obturator (BLUE)
4. Ionwand obturator (RED)

Indications

The Ultravision 5mm Trocar establishes a path of entry for instruments used in laparoscopic surgery and includes an Ionwand™ for the clearance of smoke and other particulate matter that is created during laparoscopic surgery.

Contraindications

None.

Warnings

Ultravision is not intended to be used in non-laparoscopic surgery procedures.

Use Ultravision with caution in the presence of internal or external pacemakers. Similar to electrosurgical devices, Ultravision, can affect the performance of a pacemaker. Consult the pacemaker manufacturer or hospital Cardiology Department for further information when use of Ultravision is planned in patients with cardiac pacemakers.

Similar to electrosurgical devices, Ultravision may also interfere with the performance of an internal cardiac defibrillator (ICD). If the patient has an ICD contact the ICD manufacturer for instructions before performing an electrosurgical procedure with Ultravision

The Ultravision 5mm Trocar has been designed specifically for use with the Ultravision System. Only use items supplied by Alesi Surgical Limited.

The Ultravision 5mm Trocar should only be used by surgeons with adequate training and familiarity with laparoscopic techniques.

If the pack is open or damaged, do not use and discard safely.

If there is evidence of damage to any of the components, do not use and discard safely.

Single use only – the components have not been designed for re-sterilisation and re-use of this device poses a severe risk of cross-infection and/or device malfunction.

Care should be exercised when introducing or removing sharp-edged laparoscopic instruments to lower the risk of inadvertent damage to the seal system.

Always insert the Ionwand into the patient via the Ultravision 5mm Trocar before switching on the Ultravision generator.

Insert the trocar under direct vision and avoid iatrogenic injury to abdominal structures, major arteries and veins.

As with electrosurgical instruments, contact between the Ionwand and uninsulated metallic instruments should be minimised to prevent electrostatic discharge to the operator or patient that may also transiently interfere with an ECG signal.

Repeated contact of the Ionwand with a video laparoscope may lead to temporary interference with the visual display and/or damage to the laparoscope.

Ensure that the Ultravision generator is switched off before removing the Ionwand.

Dispose of the Ultravision 5mm Trocar following the procedures for your

institution.

The seal of the Ultravision 5mm Trocar can accommodate instruments with diameter of 5mm only. Use of instruments with diameter less than 5mm can result in loss of pneumoperitoneum.

Before use

Read the Ultravision System User Manual for detailed instructions on how to set up the Ultravision System.

Insertion of the Ultravision 5mm Trocar and lonwand

1. Open the sterile pack and remove the contents using sterile field technique.
2. Remove the protective sheath from the lonwand and check the integrity (Figure 2).
3. Insert the main obturator into the trocar and ensure correct assembly. Ensure the lonwand obturator is in the trocar and luer lock is secure.
4. Make an abdominal incision that is slightly smaller than the diameter of the trocar to avoid trocar instability. Insert the trocar under direct vision.
5. Apply light and continuous, but controlled downward pressure to the obturator to introduce the trocar through the abdominal wall (Figure 3). Ensure that the entire word ULTRAVISION is visible to the laparoscope inside the peritoneum.
6. When the trocar is in the desired position, remove the main obturator and the lonwand obturator (Figure 4).
7. Fully insert the lonwand and secure in place by tightening the luer lock (Figure 5).
8. Connect the lonwand connector to the Ultravision generator (Figure 6).

During surgical procedure

9. Ensure that the lonwand is free of surgical debris and not touching patient tissue or a surgical instrument.
10. Switch on the Ultravision generator.
11. Perform the surgical procedure.
12. In the event that the proximity indicator sounds when the lonwand is not touching patient tissue or a surgical instrument, perform the following steps:
 - a. Switch off the Ultravision generator.
 - b. Remove the lonwand from the trocar.
 - c. Clean the lonwand cable and brush filaments by enclosing it within a sterile swab moistened with sterile saline and gently wiping it, moving the swab from the insulated portion of the lonwand cable towards the metal brush. Repeat until all contamination has been removed.
 - d. Re-introduce the lonwand into the trocar.
 - e. Switch on the Ultravision generator and confirm that the indicator condition has resolved.
 - f. If the alarm continues to sound, replace the lonwand.
13. Switch off Ultravision generator once surgical cutting is complete.

End of surgical procedure

14. Ensure that the Ultravision generator is off.
15. Remove the lonwand from the trocar.
16. Remove the trocar from the abdomen. If port site haemostasis is not present, appropriate technique should be used to achieve haemostasis.
17. Apply dressing to the wound site.
18. Dispose of all items safely.

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© Alesi Surgical Limited 2020
Advena Ltd., Tower Business Centre