

MT03 (12 months)

The iTind Temporarily Implanted Nitinol Device for the Treatment of Lower Urinary Tract Symptoms Secondary to Benign Prostatic Hyperplasia: A Multicenter, Randomized, Controlled Trial

Bilal Chughtai, Dean Elterman, Neal Shore, Marc Gittleman, Jay Motola, Sheldon Pike, Craig Hermann, William Terrens, Alfred Kohan, Ricardo R. Gonzalez, Aaron Katz, Jeffery Schiff, Evan Goldfischer, Ivan Grunberger, Le Mai Tu, Mark N. Alshak, and Jed Kaminetsky.

<https://pubmed.ncbi.nlm.nih.gov/33373708/>

Objective

To report the results of a multicenter, randomized, controlled trial with a temporarily implanted nitinol device (iTind; Medi-Tate Ltd, Hadera, Israel) compared to sham for the treatment of lower urinary tract symptoms secondary to benign prostatic hyperplasia.

Materials and Methods


Men 50 years or older were randomized 2:1 between iTind and sham procedure arms. A self-expanding, temporary nitinol device was placed for 5 to 7 days and an 18F Foley catheter was inserted and removed for the iTind and sham group, respectively. Patients were assessed at baseline, 1.5, 3, and 12 months post-operatively using the IPSS, peak urinary flow rate (PFR), residual urine, quality of life (QoL), and the International Index of Erectile Function (IIEF). Unblinding occurred at 3 months.

Results

A total of 175 men (mean age 61.1±6.5) participated (118 iTind vs 57 sham). 78.6% of patients in the iTind arm showed a reduction of ≥3 points in IPSS, versus 60% of patients in the control arm at 3 months. At 12 months, the iTind group reported a 9.25 decrease in IPSS ($p < 0.0001$), a 3.52ml/s increase in PFR ($p < 0.0001$) and a 1.9-point reduction in QoL ($p < 0.0001$). Adverse events were typically mild and transient, most Clavien-Dindo grade I or II, in 38.1% of patients in the iTind arm and 17.5% in the control arm. No de novo ejaculatory or erectile dysfunction occurred.

Conclusion

Treatment with the second-generation iTind provided rapid and sustained improvement in lower urinary tract symptoms for the study period while preserving sexual function.



Manufactured by Medi-Tate Ltd., 17 Hauman Street, Hadera, 3850169 Israel. Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer. Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates. | Medical devices listed may not be available for sale in all countries.