



Product Bulletin 402

New Product:

Product Update:

Obsolete Product:

Issue: Daylight Savings Time Hotfix for nStream, nStream+ and VaultStream

Date Released: 03/07/2007

Audience: Sales, Sales Support, Technical Support

Systems Affected:

ISM-020031	nStream™ SingleStream
ISM-020831	nStream™ SingleStream with DICOM
ISM-020030	nStream™ DualStream™
ISM-020830	nStream™ DualStream™ with DICOM
ISM-022031	nStream+™ SingleStream
ISM-022831	nStream+™ SingleStream with DICOM
ISM-022030	nStream+™ DualStream
ISM-022830	nStream+™ DualStream with DICOM
ISM-060050	VaultStream Server, 1500 (Pedestal)
ISM-060053	VaultStream Server, 1500 (RackMount)
ISM-060052	VaultStream Server, 2400 (Pedestal)
ISM-060054	VaultStream Server, 2400 (RackMount)

Issue:

Starting in the spring of 2007, daylight saving time (DST) start and end dates for the United States will transition to comply with the Energy Policy Act of 2005. DST dates in the United States will start three weeks earlier (2:00 A.M. on the second Sunday in March) and will end one week later (2:00 A.M. on the first Sunday in November).

Certain ISM products need to be updated so that procedures times are correct.

System Prerequisites:

- nStream V3.0, V3.1, V3.3
- nStream+ V3.4
- VaultStream V2.0, V2.1, V2.2, V2.2.1

Download:

nStream, nStream+ Microsoft HotFix KB928388:

<http://www.imagestreammedical.com/content/support/secure/tech-support/WindowsXP-KB928388-x86-ENU.zip>

VaultStream Microsoft HotFix KB928388:

<http://www.imagestreammedical.com/content/support/secure/tech-support/WindowsServer2003-KB928388-x86-ENU.zip>

Process:


Download and install the update (nStream, nStream+ process shown below):

- Unzip (if necessary) & double-click on the hot-fix executable:




Despite the Microsoft statement that you should back up your system, Image Stream Medical does not recommend this step

- Select Next.



- Select "I Agree" and click Next.



- Wait for the HotFix to update the system ...



- Click Finish.

A system reboot is not required. The process is complete.

For further questions concerning this notification please call Order Administration at (978) 486-8494.