

QuickClip™ 2 Device

Single Use Rotatable Clip Fixing Device
HX-201LR-135, HX-201UR-135, HX-201LR-135L, HX-201UR-135L, HX-201YR-135
Patient Implant Card

Patient Name: _____

Date Placed: _____

Location Placed: _____

Number of Clips Placed: _____



-This patient has received temporary mechanical endoscopic clipping device.
-Do not enter an MRI scanner room or an MR system. Failure to follow the recommended MR Conditional labeling may result in a deployed GI hemoclip dislodging from tissue or heating of tissue at the GI hemoclip location. A GI hemoclip dislodgement may result in rebleeding requiring additional intervention or surgery, serious injury, or death. If an MRI exam is needed within 6 weeks of clip placement, the patient should undergo an x-ray exam to confirm that the clips are no longer present. While most clips have passed within 6 weeks, it is possible that some clips may still be present months after the procedure.

Manufactured By
OLYMPUS MEDICAL SYSTEMS CORP.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan. Fax: (042)646-2429, Telephone: (042)642-2111

Distributed By
OLYMPUS AMERICA INC. U.S.A., Telephone: (484)896-5000
OLYMPUS CANADA INC. Canada, Telephone: (800)387-0437
OLYMPUS LATIN AMERICA, INC. U.S.A., Telephone: (305)266-2332

QuickClip™ 2 Device

Single Use Rotatable Clip Fixing Device
HX-201LR-135, HX-201UR-135, HX-201LR-135L, HX-201UR-135L, HX-201YR-135
Patient Implant Card

Patient Name: _____

Date Placed: _____

Location Placed: _____

Number of Clips Placed: _____



-This patient has received temporary mechanical endoscopic clipping device.
-Do not enter an MRI scanner room or an MR system. Failure to follow the recommended MR Conditional labeling may result in a deployed GI hemoclip dislodging from tissue or heating of tissue at the GI hemoclip location. A GI hemoclip dislodgement may result in rebleeding requiring additional intervention or surgery, serious injury, or death. If an MRI exam is needed within 6 weeks of clip placement, the patient should undergo an x-ray exam to confirm that the clips are no longer present. While most clips have passed within 6 weeks, it is possible that some clips may still be present months after the procedure.

Manufactured By
OLYMPUS MEDICAL SYSTEMS CORP.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan. Fax: (042)646-2429, Telephone: (042)642-2111

Distributed By
OLYMPUS AMERICA INC. U.S.A., Telephone: (484)896-5000
OLYMPUS CANADA INC. Canada, Telephone: (800)387-0437
OLYMPUS LATIN AMERICA, INC. U.S.A., Telephone: (305)266-2332